

# HAGMANN

*Precise • flexible • reliable*



## *Blisters and sealing lids*

*from cleanroom class ISO 7 (10'000)*



## Blister packaging solutions

Complete packaging engineering comprising:




- Optimized blister design (ergonomics and logistics)
- Appropriate punching, thermoforming and heat sealing tools
- Sealing and punching machines suitable for validation
- Standardized blister production process



## Thermoforming and cutting

Many years of experience in thermoforming

- Design and manufacture of thermoforming, sealing and cutting tools
- Production of blisters in class 7 cleanroom according to ISO 14644-1
- Cutting of sealing lids in class 7 cleanroom according to ISO 14644-1
- Printing of sealing lids
- Maximum blister sheet size 650 x 450 mm



**Hagmann Tec Ltd.**  
Bohnackerweg 6  
CH - 2545 Selzach  
Phone +41 32 641 66 41

[info@hagmann-tec.com](mailto:info@hagmann-tec.com) / [www.hagmann-tec.com](http://www.hagmann-tec.com)



Swiss Made



Reserved to technical changes /2020 / Conformity CE

