

# **HAGMANN**

*precise - flexible - reliable*




## ***Infection Control***

- Project planning
- Cleaning and disinfection devices
- Sterilizers
- Endoscope reprocessing
- Service, maintenance and validation



## Project planning

Hagmann offers comprehensive solutions for sterile processing in the AEMP:



### Project preparation

Our Sales staff record all information about premises and requirements. Regardless of the size of the facility or the project, we will take into account from the outset all the variables involved in implementing the project.

### Planning

Hagmann, as your planner, visualizes complete installation drawings taking into account the local conditions on site.

### Installation

Our service team carries out the installation, commissioning and instruction of the personnel and guarantees a professional handover from a single source.

### Training

For optimum performance of the systems, Hagmann offers customized consulting and training.

### Technical support

The Hagmann Tec AG customer service in Switzerland is ready to help you with a competent service team to solve technical problems immediately.

## Cleaning and disinfection devices

The new generation of MATACHANA MAT LD series washer-disinfectors guarantee highly efficient cleaning processes in the shortest possible time to increase productivity in the reprocessing unit for medical devices (AEMP).

➤ Maximum efficiency, minimum cycle times



# NEUE GENERATION MAT LD

INNOVATING  
HEALTHCARE

INNOVATION  
FLEXIBILITÄT  
BEDIENERFREUNDLICH  
EFFIZIENZ  
TECHNOLOGIE  
KONNEKTIVITÄT  
SICHERHEIT  
NACHHALTIGKEIT



Seamless NULLGAP installation

### ➤ HIGH LOADING CAPACITY

Up to 18 sieve baskets

### ➤ HIGH PRODUCTIVITY

Standard program 45 minutes ( LD1000 + LD500)

### ➤ LOW WATER CONSUMPTION

25 liters / cycle (LD1000)

15 liters / cycle (LD500)

### ➤ LOW SOUND PRESSURE LEVEL

< 60 dBA

### ➤ COMPACT DESIGN

Footprint 820 x 940 mm ( LD1000)

Footprint 660 x 790 mm (LD500)

## Sterilizers

The MATACHANA sterilizers were specially developed for the needs and requirements of reprocessing units for medical devices (AEMP) in hospitals, in surgical and outpatient practices and for sterilization centers of external service providers.



### Steam sterilizers

- HIGH LOADING CAPACITY**  
up to 18 STE
- FEEDING SYSTEM**  
height adjustable
- VENTURI VACUUM SYSTEM**  
Optionally with vacuum pump
- LOADING SYSTEM**  
Easy-Charge-System

### Plasma low-temperature sterilizer

**NEUE**  
**130HPO**  
**INNOVATING**  
**HEALTHCARE**

INNOVATION  
FLEXIBILITÄT  
BEDIENERFREUNDLICH  
EFFIZIENZ  
TECHNOLOGIE  
KONNEKTIVITÄT  
SICHERHEIT  
NACHHALTIGKEIT



h2o2 plasma low-temperature sterilizer

## Endoscope preparation

As the exclusive representative of Steelco for endoscope reprocessing, Hagemann Tec AG is your partner to realize a comprehensive solution.





## The company

Hagmann Tec AG is a family-owned stock corporation and was founded in 1982. Today it employs over 55 people.

Hagmann is ISO 13485 : 2016 certified. This system helps us to continuously improve and adapt to the latest standards.

## Services

For Hagmann, supporting our customers with a service is extremely important.

Our aim is to respond quickly and competently to the needs of our customers.

The comprehensive and rapid support provided by our specialists and our highly qualified and motivated staff guarantee you maximum efficiency and flawless performance of the Hagmann products .



## Our location:



**Hagmann Tec AG**  
Bohnackerweg 6 / CH-2545 Selzach Phone  
+41 32 641 66 41 / Fax +41 32 641 66 40  
info@hagmann-tec.com / www.hagmann-tec.com

